

N1081B Configurable Logic Module

Andrea Abba (abba@nuclearinstruments.eu)
CERN – Electronic Pool – 22/6/2021



N1081B / DT1081B





A pocketknife tool
in your laboratory!






The new way: N1081B

- **Multiple function** in a single NIM module
- **Replace** old function-specific NIM modules with a reconfigurable device
- **Reduce** the number of cables in your experiment
- **Simplify** the debug and verification process of the experimental setup
- **Remote reconfiguration** and data readout




22 Different functions in a single device



Remote Control & Monitor




CAEN x1081B can be configured remotely and data can be analyzed via web based interface



Remote Control & Monitor


CAEN x1081B can be configured remotely and data can be analyzed via web based interface

Section A - Ratemeter




Remote Control & Monitor

CAEN x1081B can be configured remotely and data can be analyzed via web based interface




Remote Control & Monitor


CAEN x1081B can be configured remotely and data can be analyzed via web based interface




Device I/O




Device Architecture



Section Architecture




Input discriminator




- **Bipolar discriminator** with very low selectable threshold (10mV)
- **Narrow Transition** with less than 2mV transition threshold (0 to 99%)
- **Calibrated threshold**
Non linearity $\pm 6.5\text{mV}$ $\pm 4\%$
- **Extended threshold range** from -800mV to 2.5V programmable in step of 1mV
- **Fast**
 - 130 MHz with scaler
 - 100 MHz logic function NIM
 - 80 MHz logic function NIM
 - 40 MHz synchronous functions

Input discriminator




- **Bipolar discriminator** with very low selectable threshold (<10mV)
- **Narrow Transition** with less than 2mV transition threshold (0 to 99%)
- **Calibrated threshold**
Non linearity $\pm 6.5\text{mV} \pm 4\%$
- **Extended threshold range** from -800mV to 2.5V programmable in step of 1mV
- **Fast** (up to 120 MHz with scaler, 65 MHz with logic function and 35 MHz with timer, counters, TOT)

Input gate and delay generator




- **Arbitrary shape** input signal
- **Independent** gate and delay **for each input**
- **Programmable** from 15ns to 10us, independent delay and gate
- **Resolution** +-1%+-3ns
- **Internal ring oscillator** generate pulses synchronous with the input source minimizing the input jitter
- **Temperature compensation** automatically adjust the ring oscillator frequency to generate accurate delays




Input gate and delay generator

INPUT TO
OUTPUT DELAY




- **Arbitrary shape** input signal
- **Independent** gate and delay **for each input**
- **Programmable** from 15ns to 10us, independent delay and gate
- **Resolution** $\pm 1\% \pm 3\text{ns}$
- **Internal ring oscillator** generate pulses synchronous with the input source minimizing the input jitter
- **Temperature compensation** automatically adjust the ring oscillator frequency to generate accurate delays

Output block diagram




- **TTL / NIM** input signal
- **Fast** output signal edge (<2ns)
- **Maximum rate**
100MHz NIM and 80 MHz TTL




1ns on NIM

Core functions




-  Asynchronous logic (there is no clock sampling signal in the function and the propagation delay is proportional to the number of logic element the signal is crossing)
-  Synchronous to the system 100 MHz clock (the outputs of the function are updated only on rising edge of the clock signal)

Fanout




Monitor

Activity Monitor




Logic Analyzer Monitor




Monitor

Activity Monitor



Logic Analyzer Monitor



COMPLEX PROGRAMMABLE TRIGGER TO
SEEK FOR RARE EVENTS

Trigger Settings

Mode Input

OR

Mode Output

OR

Signal Edge

FALLING

Section A

Input

1 2 3
4 5 6

ALL

Output

1 2
3 4

ALL

Section B

Input

1 2 3
4 5 6

ALL

Output

1 2
3 4

ALL

Section C

Input

1 2 3
4 5 6

ALL

Output

1 2
3 4

ALL

Section D

Input

1 2 3
4 5 6

ALL


Output

1 2
3 4


ALL

Apply


Logic functions




Scaler






SCALER (DIVISOR = 4)



Synchronous with input signal




Rate Meter




Coincidence Counter


Counts coincidences with channel 1



Counts coincidences with first arrival




Counts first arrival in respect of a specific channel



* Counts total events and channel specific coincidence

Timer

- **Up/down** time chronometer
- **Internal and external** time base (internal 100MHz)
- **Calibrated threshold** with $\pm 5\text{mV}$ $\pm 1\%$ tolerance on the full range
- **Programmable target**
- **Run and Marker (preset)** output
- 64 bit depth
- Programmable auto-reset function
- Time-gated counter



Programmable LUT

- Implement any 6 input - 4 output **logic function** for each N1081B section
- **Programmable logic table** via web interface
- Load **ASCII file** generated with any logic synthesis tool
- Web based LUT editor

Section C - Lookup Table

Start Pause File

Input Enable: 2

File Mode: Select


Select File: test2

File Name: Insert file

In 2	In 1	Out 4	Out 3	Out 2	Out 1
0	0	0	0	0	0
0	1	0	0	0	1
1	0	0	0	1	0
1	1	0	0	1	1

Page 1 of 1 View 1 - 4 of 4

Apply



Time Tagging

- **Timestamp** up to 24 detector output with
- **10ns** resolution
- Up to 50KHz long time sustainable average rate per section, up to 20 MHz peak rate.
- **Integrated discriminator** allows to direct connected detector / pre to the device
- ASCII format file downloadable from web interface
- Acquisition time limited by PC available disk space
- Up to 4 channel per section
- T0 and GATE input




		TIMESTAMP (ns)													
		CHANNEL													
1	1	8	2	5	2	7	2	4	9	2	4	8	1	2	0
2	1	8	2	5	2	7	2	5	0	2	4	8	1	2	0
3	1	8	2	5	2	7	2	5	1	2	4	8	1	2	0
4	1	8	2	5	2	7	2	5	2	2	4	8	1	2	0
5	1	8	2	5	2	7	2	5	3	2	4	8	1	2	0
6	1	8	2	5	2	7	2	5	4	2	4	8	1	2	0
7	1	8	2	5	2	7	2	5	5	2	4	8	1	2	0
8	1	8	2	5	2	7	2	5	6	2	4	8	1	2	0
9	1	8	2	5	2	7	2	5	7	2	4	8	1	2	0
10	1	8	2	5	2	7	2	5	8	2	4	8	1	2	0
11	1	8	2	5	2	7	2	5	9	2	4	8	1	2	0
12	1	8	2	5	2	7	2	6	0	2	4	8	1	2	0

Time over Threshold

- Bin size from 10ns to 1s
- Up to 1024 bin
- **Integrated discriminator** allows to direct connected detector / pre to the device
- **Online Spectrum** calculation
- ASCII format file downloadable from web interface



Time over Threshold




*measured by CAEN I-Spector device with 18x18x30mm CsI(Tl) scintillator, 10ns bit size

Time of flight

- **Bin size from 10ns to 1s**
- Up to **1024 bin**
- **Integrated discriminator** allows to direct connected detector / pre to the device
- **Online Spectrum** calculation
- ASCII format file downloadable from web interface
- Independent programmable bin size


Neutron ToF on ISIS



Time of flight


- **Bin size from 10ns to 1s**
- Up to **1024 bin**
- **Integrated discriminator** allows to direct connected detector / pre to the device
- **Online Spectrum** calculation
- ASCII format file downloadable from web interface
- Independent programmable bin size

Cadmium ToF with clearly visible resonance




Generator

- **Embedded programmable signal generator**
- **Periodic and Poisson distributed events**
- 1Hz to 20MHz generation
- **10ns to 1s pulse width**
- **Low / High constant generator**




Section bypass



- **Daesy chain** section to extend the number of input of each function





Propagation delay






External clock




Localization in crowded crate






User Manual



Timing diagram of every function is every possible configuration



Online help



Configuration Memory

SAVE

File name

Config_

SAVE

✕

🏠

LOAD

Config_test

Config_Demo

Config_rma_pulse

LOAD

✕

🏠

A hand is raised in the air, palm facing forward, against a blurred background of a crowd. The hand is in sharp focus, while the background is out of focus, showing other people and what appears to be a stage or presentation area.

Thank you for your attention

Any question/curiosity?